

Promiseb[®] Topical Cream

For Topical Dermatological Use Only

Rx only

For External Use Only

Product Description

Promiseb[®] Topical Cream is an off-white, steroid-free, fragrance-free, water-based emulsion.

Indications for Use:

Under the supervision of a healthcare professional, Promiseb Topical Cream is indicated to manage and relieve the signs and symptoms of seborrhea and seborrheic dermatitis such as itching, erythema, scaling and pain. Promiseb Topical Cream helps to relieve dry waxy skin by maintaining a moist wound & skin environment, which is beneficial to the healing process.

Directions for Use:

Apply Promiseb Topical Cream to the affected skin areas 2 to 3 times per day (or as needed), and massage gently into the skin. If the skin is broken, cover Promiseb Topical Cream with a dressing of choice.

Ingredients:

Promiseb Topical Cream is comprised of Purified Water, Isohexadecane, Butyrospermum parkii, Pentylene glycol, Ethylhexyl palmitate, Cera alba, PEG-30 Dipolyhydroxystearate, Bisabolol, Polyglyceryl-6 polyricinoleate, Tocopheryl acetate, Hydrogenated castor oil, Acifructol complex, Butylene glycol, Magnesium sulfate, Piroctone olamine, Allantoin, Magnesium stearate, Disodium EDTA, Vitis vinifera, Ascorbyl tetraisopalmitate, Glycyrrhetic acid, Propyl gallate, and Telmestaine.

Caution:

The use of Promiseb Topical Cream is contraindicated in any patient with a known history of hypersensitivity to any of the ingredients. Promiseb Topical Cream does not contain milk, wheat, peanut or animal derivatives. Promiseb Topical Cream does contain shea butter (Butyrospermum parkii), a derivative of shea nut oil (not peanut oil). Patients with a known allergy to nuts or nut oils should consult their physician before using this topical preparation.

How Supplied

30g tube, 67857-803-30

60g (2 x 30g Tubes), 67857-803-60

To Open: Puncture seal with pointed end of cap.

Important: The opening of this product is covered by a metal seal. Do not use if seal has been punctured or is not visible.

Store at controlled room temperature 68° to 77°F (20° to 25°C), excursions permitted between 59° and 86°F (15° and 30°C).

Distributed by: Promius Pharma, LLC, Princeton, NJ 08540

PROMIUS[®]
P H A R M A

Made in Italy

Federal law restricts this device to sale by or on the order of a physician or properly licensed practitioner.